TRUEtest™ Blood Glucose Test Strips Instructions for Use

Intended Use

TRUEtest™ Blood Glucose Test Strips are used only with TRUEresult™ and TRUE2go™ Meters to quantitatively measure whole blood glucose. TRUEtest Strips can be used for self-testing at home or for professional use. TRUEtest is for *in vitro* diagnostic use only.

Test Principle

The TRUEtest Strip is a plastic strip containing chemicals and electrodes. When inserted into a TRUEresult or TRUE2go Meter, glucose is measured using amperometric technology employing a glucose dehydrogenase-PQQ reaction. When whole blood or Control is drawn into the Sample Tip of the Strip, glucose in the sample reacts with the chemicals and produces an electrical current. The Meter measures the current and calculates the amount of glucose. The result is displayed as a plasma value.

Chemical Composition

Glucose dehydrogenase-PQQ (*Acinetobacter calcoaceticus*), mediators, buffers and stabilizers.

① **Contact End** End of Strip inserted into Meter.

② **Sample Tip** Pointed end of Strip where sample is drawn into Strip.

Top of Test Strip ②

①

Correct Incorrect

Caring for Strips

- Strips must be kept in original vial with vial cap tightly sealed. NEVER transfer Strips from one vial to another.
- Use Strip quickly after removing from vial. Recap vial right away. Strips left outside of vial too long give an error message.
- Write date opened on Strip vial label when removing the first Strip. Discard all unused Strips in vial after date printed on the Strip vial label, or 4 months after date opened, if either date has passed.
- Store Strip vial in a dry place at room temperature below 86°F. DO NOT REFRIGERATE OR FREEZE. Do not store in bathroom or kitchen. Do not expose to extreme heat or cold, direct sunlight or high humidity for any length of time.
- Do not bend, cut, or alter Strips in any way.

Important Information

- Strips are for *in vitro* testing only. Do not consume.
- Use TRUEtest Strips only with TRUEresult and TRUE2go Meters and TRUEtest Glucose Control. Using other meters or Controls may give inaccurate results.
- Do not change treatment plan based on the results from the TRUEtest Strips and TRUEresult/TRUE2go Meters without the advice of a Doctor or Healthcare Professional.
- Discard used Strips and lancets into an appropriate container.
- Using Strips past written or printed dates on the Strip vial label may cause inaccurate results.
- Discard any Strips or vials that appear damaged.
- NEVER use serum, plasma, or clotted blood for testing. Use fresh, capillary whole blood or venous whole blood drawn into an EDTA (purple top) or heparin (green top) vacutainer tube. Mix well before sampling.
- When using the forearm for blood sample:¹
 - Check with your Doctor or Healthcare Professional to see if forearm testing is right for you.
 - Results from forearm are not always the same as results from finger.
- Use finger instead of forearm for more accurate results:
 - Within 2 hours of eating, exercise, or taking insulin,
 - If your blood sugar may be rising or falling rapidly or your routine results are often fluctuating,
 - If you are ill or under stress,
 - If your forearm test results do not match how you feel,
 - If your blood sugar may be low or high,
 - If you do not notice symptoms when blood sugar is low or high.

Quality Control Testing

There are two quality control tests to let you know that the System is working properly. An automatic self-test is performed each time a TRUEtest Strip is inserted into a TRUEresult or TRUE2go Meter. Upon inserting a Strip into the Test Port, if all segments appear and the Apply Sample Symbols appear in the Display, the Meter is working properly and can be used for testing.

Important Information: There are three levels of TRUEtest Glucose Control available that contain known amounts of glucose. It is important to perform Control Tests with more than one level of Control to assure your System is working properly and your testing technique is good. For more information on obtaining different levels of Control, call 1-800-803-6025 or 1-954-677-4599. See TRUEtest Glucose Control Instructions for Use or TRUEresult/TRUE2go Owner's Booklet for more information on Quality Control testing.

Blood Glucose Testing

- Note to Healthcare Professionals: Venous whole blood collected in EDTA (purple top) or heparin (green top) vacutainer tubes may be used for testing. Mix blood thoroughly before testing.
1. Wash area to be lanced, dry.
 2. Allow Meter and Strips to sit at room temperature for 30 minutes. If opening vial for the first time, write date opened on vial label.
 3. Check written date and printed date on Strip vial label. Do not use if either date has passed. Discard vial and test with new vial.
 4. Remove one Strip from vial. Recap vial right away.
 5. Insert Contact End of Strip into Test Port of Meter. Meter turns on. Do not remove Strip from Meter until testing is finished.
 6. Obtain blood drop.
 7. After Apply Sample Symbol appears in Display, place Sample Tip of Strip gently against blood drop and hold. Allow Strip to draw in blood (0.5 µL) until dashes appear in Display. Remove Strip from drop. Dashes appear as Meter is testing. If dashes do not appear in Display within 5 seconds after touching Sample Tip to drop, see Troubleshooting in the TRUEresult / TRUE2go Owner's Booklet.
 8. Result is displayed. Record result.
 9. Hold Meter with Strip pointing down. Press Strip Release button to discard Strip into appropriate container.

Expected Results for people without diabetes¹:

	Plasma Blood Glucose Result
Before eating	<100 mg/dL
Two hours after meals	<140 mg/dL

A Doctor or Diabetes Healthcare Professional determines personal target glucose ranges. If you are having symptoms that suggest your glucose is too low or too high, contact your Doctor or Diabetes Healthcare Professional right away. If comparing results using TRUEtest strips to laboratory test results, perform a fingerstick blood test within 15 minutes of the laboratory test. Results are considered accurate when they are within 20% of laboratory results. If you have eaten recently, results using TRUEtest Strips can be up to 70 mg/dL higher than laboratory results.²

Troubleshooting

- If your result is unusually high or low or doesn't match the way you feel, perform a Control Test (see **Quality Control Testing**).
- If the Control Test is within range:
- Read **Blood Glucose Testing** again.
 - Recheck your results with a new TRUEtest Strip.
- If the results still do not match the way you feel, check with your Doctor or Healthcare Professional before changing your treatment program.
- If the results are not within range:
- Check the Expiration Dates. Do not use if past either written date or date printed on Test Strip vial or Control bottle. Test with new Strips/Control.
 - Check for error messages. If an error message appears, follow the Actions in the Display Message Section of the Owner's Booklet
 - Check your testing technique. Perform another Control Test.
 - Check the temperature. Allow System to reach room temperature 59-86°F before testing.

Limitations³

WARNING!

- **Peritoneal dialysis patients receiving dialysis solutions containing Icodextrin (e.g., Extraneal®, Icodial®) should not use the TRUEresult or TRUE2go System. The dialysis solution may falsely raise glucose results.**
 - **Injection or infusion of solutions containing galactose or maltose (includes some human immunoglobulin preparations) may falsely raise glucose results.**
 - **Blood concentrations of galactose > 10 mg/dL, maltose > 12.5 mg/dL, maltotriose > 20 mg/dL and maltotetraose > 10 mg/dL may falsely raise glucose results.**
 - **Do not use TRUEresult or TRUE2go Systems during a xylose absorption test. This may falsely raise glucose results.**
 - Blood samples containing high uric acid concentration (≥ 9 mg/dL) at glucose levels < 240 mg/dL may be detected as control samples by the TRUEresult Meter, but not the TRUE2go Meter.
- The following will not affect accurate results:
- Testing at altitudes up to and including 10,150 feet.
 - Hematocrit levels between 20% and 60%.
- The following may affect accurate results:^{4,5,6}
- Hyperglycemia with hyperosmolarity, with or without ketosis.
 - Thickened blood caused by dehydration.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock.⁷ Inaccurate results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.

FOR CONSUMERS • Performance Characteristics⁴

Accuracy: Diabetes experts have suggested that glucose meters should agree within 15 mg/dL of a laboratory method when the glucose concentration is less than 75 mg/dL, and within 20% of a laboratory method when the glucose concentration is 75 mg/dL or higher. TRUEresult and TRUE2go were tested by users at diabetes clinics, large urban hospitals, and diabetes care centers. The table below shows how often user TRUEresult and TRUE2go fingertip values can achieve these goals. The fingertip data were compared to parallel results obtained on a Yellow Springs Instrument (YSI) Model 2300.

TRUEresult Finger Sample - <75 mg/dL (user finger vs. YSI)			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
16/25 = 64%	23/25 = 92%	25/25 = 100%	
TRUEresult Finger Sample - ≥ 75 mg/dL (user finger vs. YSI)			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
169/321 = 53%	271/321 = 84%	310/321 = 97%	319/321 = 99%
TRUE2go Finger Sample - <75 mg/dL (user finger vs. YSI)			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
15/22 = 68%	22/22 = 100%	22/22 = 100%	
TRUE2go Finger Sample - ≥ 75 mg/dL (user finger vs. YSI)			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
171/310 = 55%	255/310 = 82%	299/310 = 96%	309/310 = 99.7%

The table below shows how often user TRUEresult and TRUE2go forearm values achieve these goals when users' glucose values are not fluctuating.

TRUEresult Forearm Sample - <75 mg/dL (user forearm vs. user finger)			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
17/25 = 68%	23/25 = 92%	25/25 = 100%	
TRUEresult Forearm Sample - ≥ 75 mg/dL (user forearm vs. user finger)			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
193/319 = 61%	285/319 = 89%	311/319 = 97%	314/319 = 98%
TRUE2go Forearm Sample - <75 mg/dL (user forearm vs. user finger)			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
14/19 = 74%	19/19 = 100%	19/19 = 100%	
TRUE2go Forearm Sample - ≥ 75 mg/dL (user forearm vs. user finger)			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
203/311 = 65%	269/311 = 86%	299/311 = 96%	306/311 = 98%

FOR HEALTHCARE PROFESSIONALS⁴

Performance Characteristics – TRUEresult System

Accuracy: TRUEresult accuracy was assessed against the Yellow Springs Instrument (YSI) Model 2300. Studies were conducted at 4 clinical sites by Healthcare Professionals.

Fingertip Capillary Blood (mg/dL) (ISO 15197:2003 data finger vs. YSI):⁴ N=314

Slope	y-intercept	r ²	Range
1.01	-2.19	0.98	24-549

Venous Blood (mg/dL):⁴ N=236

Slope	y-intercept	r ²	Range
1.11	-4.40	0.99	34-463

99.4% of Healthcare Professional (HCP) TRUEresult fingertip values fell within 15 mg/dL of the YSI results at glucose levels <75 mg/dL and within 20% at glucose levels ≥ 75 mg/dL.

Fingertip Capillary Blood (HCP finger vs. YSI)

<75 mg/dL			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
28/44 = 64%	42/44 = 95%	44/44 = 100%	
≥ 75 mg/dL			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
139/270 = 51%	227/270 = 84%	260/270 = 96%	268/270 = 99%

98.8% of Healthcare Professional (HCP) TRUEresult forearm values fell within 15 mg/dL of the fingertip results at glucose levels <75 mg/dL and within 20% at glucose levels ≥ 75 mg/dL when users' glucose values are not fluctuating.

Forearm Capillary Blood (HCP forearm vs. HCP finger)

<75 mg/dL			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
12/26 = 46%	22/26 = 85%	25/26 = 96%	
≥ 75 mg/dL			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
196/318 = 62%	286/318 = 90%	311/318 = 98%	315/318 = 99%

Precision: Precision describes the variation between results. Precision results were performed in a laboratory.

Blood (Within Lot):⁴ N=100

mg/dL	44	77	130	214	319
CV%	2.5	3.5	3.2	4.6	3.9

Blood (Within Vial):⁴ N=10

mg/dL	44	78	130	217	319
CV%	2.0	2.2	2.1	3.2	2.4

Glucose Control:⁴ N=100

mg/dL	45	104	301
CV%	1.7	1.6	4.1

Performance Characteristics – TRUE2go System

Accuracy: The accuracy of TRUE2go was assessed against the Yellow Springs Instrument (YSI) Model 2300. Studies were conducted at 4 clinical sites by healthcare professionals.

Fingertip Capillary Blood (mg/dL) (ISO 15197:2003 data finger vs. YSI):⁴ N=332

Slope	y-intercept	r ²	Range
0.95	2.58	0.98	45-498

Venous Blood (mg/dL):⁴ N=236

Slope	y-intercept	r ²	Range
1.11	-3.52	0.99	34-463

100% of Healthcare Professional (HCP) TRUE2go Meter fingertip values fell within 15 mg/dL of the YSI results at glucose levels <75 mg/dL and within 20% at glucose levels ≥ 75 mg/dL.

Fingertip Capillary Blood (HCP finger vs. YSI)

<75 mg/dL			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
30/50 = 60%	48/50 = 96%	50/50 = 100%	
≥ 75 mg/dL			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
152/289 = 53%	254/289 = 88%	278/289 = 96%	289/289 = 100%

99.4% of Healthcare Professional (HCP) TRUE2go Meter forearm values fell within 15 mg/dL of the Healthcare Professional finger results at glucose levels <75 mg/dL and within 20% at glucose levels ≥ 75 mg/dL.

Forearm Capillary Blood (HCP forearm vs. HCP finger)

<75 mg/dL			
± 5 mg/dL	± 10 mg/dL	± 15 mg/dL	
18/19 = 95%	19/19 = 100%	19/19 = 100%	
≥ 75 mg/dL			
$\pm 5\%$	$\pm 10\%$	$\pm 15\%$	$\pm 20\%$
202/311 = 65%	270/311 = 87%	300/311 = 96%	308/311 = 99%

Precision: Precision describes the variation between results. Precision results were performed in a laboratory.

Blood (Within Lot):⁴ N=100

mg/dL	47	77	138	211	312
CV%	3.0	2.7	3.1	4.7	4.4

Blood (Within Vial):⁴ N=10

mg/dL	48	77	139	209	312
CV%	1.7	1.9	2.1	3.1	3.5

Glucose Control:⁴ N=100

mg/dL	45	102	298
CV%	1.7	1.7	2.8

References/Referencias

1. American Diabetes Association Standards of Medical Care in Diabetes, *Diabetes Care*, January 2008, pgs. S12-S54.
2. Joslin Diabetes Center. *Goals for Blood Glucose Control*. [Electronic Version]. Retrieved June 2, 2008: www.joslin.org/Beginners_guide_523.asp.
3. Larson-Cann U: Differences between capillary and venous blood glucose during oral glucose tolerance tests. *Scand J Clin Lab Invest* 36:805-808, 1976.
4. Data on file.
5. Young, D.S. Effects of Drugs of Clinical Laboratory Tests, Third Edition. Washington, D.C.: AACCC Press (1990), pp. 3-168, 3-180.
6. Sacks, D.B., "Carbohydrates," Burtis, C.A., and Ashwood, E.R. (ed.), *Tietz Textbook of Clinical Chemistry*. Philadelphia: W.B. Saunders Company (1994), 959.
7. Atkins, S.H.; Dasamhapatra, A.; Jaker, M.A.; Chorost, M.J.; Redd, S., *Fingerstick Glucose Determination in Shock*. *Annals of Internal Medicine*, 114:1020-1024 (1991).
8. International Organization for Standardization. *In vitro diagnostic test systems. Requirements for blood-glucose monitoring system for self-testing in managing diabetes mellitus*. Reference number ISO 15197:2003 (E). Geneva: International Organization for Standardization, 2003.

Additional Information: See the Owner's Booklets for more detailed instructions. Call Home Diagnostics, Inc. at 1-800-803-6025 (USA) or 1-954-677-4599 for assistance. For medical assistance, call your Doctor or Diabetes Healthcare Professional.

Manufactured by:
HOMEdiagnostics™
 Fort Lauderdale, FL 33309 U.S.A.
 © 08/2009 Home Diagnostics, Inc. TRUEtest, TRUEresult, TRUE2go and the Home Diagnostics logo are trademarks of Home Diagnostics Inc. E3HD103 Rev. 6

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Referencias

- American Diabetes Association Standards of Medical Care in Diabetes. *Diabetes Care*, January 2008, págs. S12-S54.
- Joslin Diabetes Center. *Goals for Blood Glucose Control*. [Versión electrónica]. Copiado el 02 de junio de 2008: www.joslin.org/Beginners_guide_523.asp.
- Larson-Gunn U. Differences between capillary and venous blood glucose during oral glucose tolerance tests. *Scand J Clin Lab Invest* 36:805-808, 1976.
- Información archivada.
- Young, D.S. Effects of Drugs of Clinical Laboratory Tests, tercera edición. Washington, D.C.: AACC Press (1990), págs. 3-168, 3-180.
- Sacks, D.B., "Carbohydrates", Burris, C.A. y Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry, Philadelphia: W.B. Saunders Company (1994), 959.
- Atkins, S.H.; Alkns, S.H.; Dasmanahapatra, A.; Jaker, M.A.; Charosi, M.L.; Radd, S., Fingerstick Glucose Determination in Shock. *Annals of Internal Medicine*, 114:1020-1024 (1991).
- Young, D.S. Effects of Drugs of Clinical Laboratory Tests, tercera edición. Washington, D.C.: AACC Press (1990), págs. 3-168, 3-180.
- International Organization for Standardization. *In vitro diagnostic test systems. Requirements for blood-glucose monitoring system for self-testing in managing diabetes mellitus*. Reference number ISO 15197:2003 (E). Geneva: International Organization for Standardization, 2003.

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