

TRUEread[®] Blood Glucose Monitoring System

What is the TRUEread[®] Blood Glucose Testing System?

Used to test blood glucose levels, the TRUEread[™] meter is easy to handle and simple to use; our TRUEread[®] meter is perfect for patients on Medicare and Medicaid plans.

What strips do I use for my TRUEread[®] meter?

TRUEread[™] meters use [TRUEread[™] Blood Glucose Test Strips](#).

How do I set up my TRUEread[®] meter?

See [instructional video](#) and refer to your [owner's booklet](#) included with your meter.

How do I calibrate my TRUEread[®] meter?

See [instructional video](#) and refer to your [owner's booklet](#) included with your meter.

How do I use the TRUEread[®] System?

See [instructional video](#) and refer to your [owner's booklet](#) included with your meter.

Can the TRUEread[®] System be used for alternative site testing?

Yes. TRUEread[®] is suitable for fingertip and forearm testing.

Are special lancets required when testing with the TRUEread[®] meter?

No. We recommend the TRUEdraw[™] lancing device for greater comfort and control to be used with our single-use sterile lancets.

I am on peritoneal dialysis, can I use this system?

Yes. The TRUEread[®] System is a great choice for you. Please consult your test strip instructions for use for limitations.

Does the TRUEread[®] meter have memory? Does the meter store my results?

Yes. TRUEread[®] has a 200-test memory with uploading capability.

Can I upload my meter data into a data management system?

Yes. The TRUEread[®] meter uses [TRUEmanager[™] Diabetes Management Software](#). Follow the link to download the software or to request the software on CD. A SmartData Cable may be purchased separately.

Is the TRUEread[®] meter covered by a lifetime warranty?

Yes. We are committed to ensuring that your meter and testing requirements are fulfilled. Our meters are covered by a lifetime warranty.

What kind of batteries does the TRUEread[®] meter use?

The TRUEread[®] meter uses a non rechargeable 3V lithium battery (#CR2032).

Where can I find the TRUEread[®] System?

Our TRUE portfolio of products are available at your favorite retail store or mail service provider. [Where To Buy](#)